





NEWS

Endoscopic wire trays...

...for transportation, storage, cleaning and sterilization!



Advantages

Safety

- safety storage for telescopes during cleaning or sterilization.
- external contours are polished and sharp-free.


Simple using

- Interchangeable silicone parts.
- Special lock to avoid self-opening.
- no damage of packaging material.

Consisting of

- a lot of tray types.
- fine wire netting for the cleaning process.

Dimensions BxTxH	Article N°
Tray for 1 telescope	
290 x 80 x 52 mm	64-21400-01
460 x 80 x 52 mm	64-21400-02
670 x 80 x 52 mm	64-21400-03
Tray for 1 telescope with fixture for fiber optic guide including basket for small parts	
290 x 80 x 52 mm	64-21400-04



Design

- chrome-nickel steel 18/10, material 1.4301
- electrolytic polished
- sharp-free finished
- closed wire loop
- detachable lid with special lock
- frame: Ø 5 mm.
- bottom and lateral wire: Ø 2 mm.