

STERIKING®



Rotosealer RS-220 A

- ▶ All-day ready for use, simple operation
- ▶ Safety tested and approved by EN, VDE
- ▶ Microprocessor-controlled for reliable function
- ▶ Symbolic language according to EN 980
- ▶ CE-mark function
- ▶ Different printing & counter options
- ▶ Multi-language programming
- ▶ Daily validation with the seal check function
- ▶ Automatic maintenance control
- ▶ Shock watch system for delivery safety


WIPAK MEDICAL

Rotosealer RS-220 A

Simultaneous Sealing & Printing of Sterilization Packages

Continuous working unit for quick and reliable closing and marking of all heat-sealable sterilization packages.

Optimum sealing reliability

- Motor start-up protection in the event of 5 °C deviation from the set-point temperature ensures correct sealing.
- Automatic on/off switch function.
- Automatic energy-saving standby function.
- The compression rollers ensure tight closure of packages with or without gussets.
- Adjustable in-feed guiding plate provides predefinition of the seal distance (5 - 30 mm) from the package edge.



Daily validation with the Seal Check function: temperature, pressure and sealing time can be printed and recorded.

Technical data	
Seal width	12 mm (2 x 5 mm)
Sealing speed	10 m/min
Sealing temperature	80 - 220 °C
Temperature tolerance	±2 %
Power consumption	500 W
Power supply	230/115 V 50/60 Hz
Weight	23 kg
Dimensions	620 x 260 x 250 mm
Separate worktop	625 x 225 x 60 mm

Warranty

- 24 months from the date of purchase.

Shock watch label

- Enables recognition of transport damage.
- Adhered to the transport carton.

Detailed Operating Manual and List of Spare Parts come with each unit.



Printer function

- The desired data can be easily fed in and modified using the keyboard
- Sterilization and expiry dates updated automatically
- Alphanumerical lot code
- User identity number
- Symbols according to EN 980
- Freely selectable text and catalogue memory
- Programmable font style and character spacing
- CE-mark with alphanumerical characters up to 20 digits
- Adjustable print starting point
- Additional text up to 20 digits
- Print rotatable through 180°
- Printer can be switched off
- The standard ink ribbon cassette is easy to replace on the front side of the unit

Integrated counter

- Automatic batch and piece counter
- Operating hours counter
- Easy and simple to read on the display and easy to reset using the keyboard

Control panel

- Digital display easy to read
- Simple input and modification of data
- Keyboard provides clear push buttons with symbols
- Multi-language programming selection: English, German, French, Spanish and Italian
- Display shows the set-point and actual temperatures
- Clock and calendar function
- The input data is preserved after the unit is switched off



We care that you pack safely!