

Chemical indicator strips

The Steriking® chemical indicator strips for monitoring process conditions in each sterilization cycle. All indicators are corresponding to critical parameters in the sterilization cycle

STERIKING® AC250 STEAM Strips

AC	STERIKING® AC250 121°C, 15 min; 134°C, 3.5 min ISO 11140-1:2005, Class 4	Contents _____
	STEAM processed = "A" and ● turn black / schwarz / noir / nero / negro / черный	Date _____
		Initials _____

AC	STERIKING® AC250 121°C, 15 min; 134°C, 3.5 min ISO 11140-1:2005, Class 4	Contents _____
	STEAM processed = "A" and ● turn black / schwarz / noir / nero / negro / черный	Date _____
		Initials _____

After sterilization

STERIKING® AC250/500 Double STEAM Strips

STERIKING® AC DOUBLE STRIP 121°C, 15 min; 134°C, 3.5 min	● turn black / schwarz / noir / nero / negro / черный	STERIKING® AC DOUBLE STRIP 121°C, 15 min; 134°C, 3.5 min
ISO 11140-1:2005, Class 4 STEAM processed = black / schwarz / noir / nero / negro / черный		STERIKING® AC DOUBLE STRIP 121°C, 15 min; 134°C, 3.5 min

STERIKING® AC DOUBLE STRIP 121°C, 15 min; 134°C, 3.5 min	● turn black / schwarz / noir / nero / negro / черный
ISO 11140-1:2005, Class 4 STEAM processed = black / schwarz / noir / nero / negro / черный	

After sterilization

STERIKING® EO250 EO GAS Strips

EO	STERIKING® EO250 EO 600 mg/l, 37°C, 60% RH, 100 min	Contents _____
	ISO 11140-1:2005, Class 5 EO processed = green / grün / vert / verde / ЗЕЛЕНЬИЙ	Date _____
		Initials _____

EO	STERIKING® EO250 EO 600 mg/l, 37°C, 60% RH, 100 min	Contents _____
	ISO 11140-1:2005, Class 5 EO processed = green / grün / vert / verde / ЗЕЛЕНЬИЙ	Date _____
		Initials _____

After sterilization

STERIKING® FO250/500 Double FORMALDEHYDE Strips

FORM	STERIKING® FO DOUBLE STRIP ISO 11140-1, Class 1	processed = green / grün / vert / verde / ЗЕЛЕНЬИЙ
	FORM processed = green / grün / vert / verde / ЗЕЛЕНЬИЙ	STERIKING® FO DOUBLE STRIP ISO 11140-1, Class 1

STERIKING® FO DOUBLE STRIP ISO 11140-1, Class 1	processed = green / grün / vert / verde / ЗЕЛЕНЬИЙ
FORM processed = green / grün / vert / verde / ЗЕЛЕНЬИЙ	

After sterilization

STERIKING® Gas Plasma Strips

STERIKING® GAS PLASMA STRIP	REF PL500 LOT1945946
Charge: H2O2 PLASMA	Name: _____
Date: processed : blue/blau/bleu/blue/azul/СИНИЙ	Other: _____

STERIKING® GAS PLASMA STRIP	REF PL500 LOT1945946
Charge: H2O2 PLASMA	Name: _____
Date: processed : blue/blau/bleu/blue/azul/СИНИЙ	Other: _____

After sterilization